











REF P9700B1HW 0,30 kg - 0,65 lbs		REF P9750F1HW 0,23 kg - 0,50 lbs		REF P9700B1HW 0,30 kg - 0,65 lbs		REF P9750F1HW 0,23 kg - 0,50 lbs	
REF P9701B1HW 0,60 kg - 1,3 lbs		REF P9751F1HW 0,60 kg - 1,3 lbs		REF P9702B1HW 0,60 kg - 1,3 lbs		REF P9752F1HW 0,68 kg - 1,5 lbs	
REF P9703B1HW 0,45 kg - 1 lbs		REF P9753F1HW 0,42 kg - 0,9 lbs		REF P9704B1HW 0,45 kg - 1 lbs		REF P9754F1HW 0,53 kg - 1,1 lbs	
REF P9705B1HW 2,1 kg - 4,5 lbs		REF P9755F1HW 1,8 kg - 3,9 lbs		REF P9706B1HW 1,80 kg - 4 lbs		REF P9757F1HW 2,10 kg - 4,6 lbs	
REF P9708B1HW 0,95 kg - 1,95 lbs		REF P9759F1HW 0,79 kg - 1,74 lbs		REF P9709B1HW 1,25 kg - 2,75 lbs		REF P9760F1HW 0,90 kg - 1,98 lbs	
REF P9710B1HW 0,31 kg - 0,68 lbs		REF P9762F1HW 0,30 kg - 0,66 lbs		REF P9711B1HW 0,33 kg - 0,72 lbs		REF P9716F1HW 0,26 kg - 0,57 lbs	
REF P9712B1HW 0,31 kg - 0,68 lbs		REF P9763F1HW 0,75 kg - 1,65 lbs		REF P9713B1HW 0,33 kg - 0,72 lbs		REF P9716F1HW 0,26 kg - 0,57 lbs	
REF P9714B1HW 0,13 kg - 0,28 lbs		REF P9764F1HW 0,10 kg - 0,22 lbs		REF P9715B1HW 1,25 kg - 2,75 lbs		REF P9716F1HW 0,26 kg - 0,57 lbs	

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EN

Medical device of class I according MDR (EU) 2017/45.
These devices are intended to be handled by anyone able to read and understand the instructions for use.

INTENDED USE OF PRODUCT

A range of multipurpose appliances that can be used both for positioning and for prevention and treatment help of bedsores.

INDICATION

Their universal nature enables them to be used in various positions, supporting the lower limbs in the semi-fowler position, or with the patient on the side or propping up the legs. The choice of device must be approved by the healthcare professional.

CONDITIONS OF USE

Devices with an integral POLYMAILLE® / POLYMAILLE-HD cover can be used for single or multiple patients. They can be reconditioned in accordance with the product instructions and the cleaning and disinfection instructions.

The first preventive measure against bedsores: the body should be raised as frequently as possible to limit compression of the blood vessels, thereby enabling a correct circulation of the blood.

CONTRA-INDICATIONS

Do not use the device beyond the intended purposes use. Do not use if you are allergic to polyurethane. For any problem, consult your healthcare provider.

SIDE EFFECTS

There are no undesirable side effects. However, the patient may experience the effects inherent to immersion: possible (limited) reduction in perception and mobility. Increased sensation of heat.

DEVICE PERFORMANCE

Enable to relieve support areas through immersion or discharge. Helps maintaining positioning in bed. Improves tissue oxygenation in the contact zone.

INFORMATION FOR HEALTHCARE PROFESSIONAL

Use positioning aids or pressure relief devices in patient with bedsores. The choice of device must be approved suitable for the patient's environment. Ensure continuity of treatment.

STORAGE, HANDLING, TRANSPORT

- Preferably flat. Do not expose close to a significant heat source.

- Imperatively in a dry and temperate place protected from heat and light in connection bag.

- Handle the device with care to avoid damaging the cover.

- Transport in a transport bag or in its box.

- For any handling / transport, please follow these instructions.

PRECAUTIONS OF USE

Do not use without cover body between the patient and the device.

- Be careful not to damage the cover. Frequent degradations observed: cigarette burns, pet scratches, needle holes, cuts with scalpels or sharp objects, damage caused by the use of equipments with sharp edges.

- We recommend to provide an additional protective cover.

- Do not use spare covers, accessories or components other than those supplied by SYSTAM®.

- Refer to 'Cleaning & Disinfection' instructions.

- If you observe deterioration of the product (soam deformation or leakage) or bottoming effect, contact your dealer.

- Report to your caregiver as soon as possible any abnormal event such as pain or redness or whitening of support points.

- Any serious incident that has occurred in connection with the device should be notified to the manufacturer and to the local competent authority.

TECHNICAL SPECIFICATIONS

Microbeads: cushion material: expanded polyurethane (PU).

Fibers cushion material: polyester (PES) fibers.

Inside case material: 100% cotton (CO).

COVERING LIFETIME

The lifespan of these products is 5 years when used in

strict accordance with the care and maintenance instructions. The service life can be longer if the product is maintained with care and considerably shortened by extreme or improper use. This does not constitute an additional guarantee on these products.

GUARANTEE

Two years against any manufacturing defaults, from the date of purchase by the customer. Small imperfections in the manufacturing process in no way affect the reason for effectiveness of the product and cannot be a reason for replacement of the product or for a return of the device or cover. In the event of default, please contact your provider, and give them the product lot number.

SAFETY & ENVIRONMENT

Do not burn or dump. Respect the recycling channels in place in your country. Do not ingest. 0% PVC foam. Comply with fire standards EN 597-1 & 597-2.

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